

Autoreduction setup forms

Mathieu Doucet
Oak Ridge National Laboratory

ORNL is managed by UT-Battelle, LLC for the US Department of Energy



U.S. DEPARTMENT OF
ENERGY

Plan for the next few weeks


Test environment:

New RHEL8 machines are being set up so we can install them together

Topics to cover:

1. General overview
2. Workflow manager and DASMON listener – Installation & maintenance [this presentation]
3. Web monitor – Installation and maintenance
4. Autoreduction service – Installation and maintenance
5. **Autoreduction setup through webmon**
6. The IHC call – when things go wrong & recovery strategies
7. Vision for the future – what I would do differently

Post-Processing Architecture



Setup accessible through web monitor

The code is in data_workflow/reporting/reduction

User submits form

- Most recent parameter values are stored in the DB
- Form allows user to see/modify parameters

Process form

- Store new values in DB
- Compose json payload
- Submit request to AR through AMQ



Write script

- AR service receives request
- Template is used to generate script
- /SNS/CNCS/shared/autoreduce/reduce_CNCS.py.template

CNCS Run 346752

home > cncs > ipts-24051 > run 346752 live monitoring: status | runs | PVs
previous | next

Run title Sample=1, ITEMS=69411, Ei=3.32, T=400
Run start July 14, 2020, 8:50 a.m.
Run end July 14, 2020, 11:16 a.m.
Duration 8743.71679688
Total counts 18523892
Proton charge 1.19999897562e+13



Submit for post-processing: [catalog](#) | [reduction](#) | [all post-processing](#) | [setup](#)

The setup form

<https://monitor.sns.gov/reduction/cncs/>

- Code is in `data_workflow/reporting/reduction/forms.py`
- Inherit from a common base class `BaseReductionConfigurationForm`
- Values are stored as `ReductionProperty` model objects
- New params can be added to the form and will be created in the DB automatically
- An html template also needs to be modified:
 - `data_workflow/reporting/templates/reduction/configuration_cncs.html`

<input type="checkbox"/>	22	cncs	grouping	powder
<input type="checkbox"/>	21	corelli	useCC	True
<input type="checkbox"/>	20	corelli	plot_requests	{('Minimum': '-0.1', 'PerpendicularTo': '[0,K,0]', 'Maximum': '0.1')}
<input type="checkbox"/>	19	corelli	vanadium_flux_file	/SNS/CORELLI/IPTS-21454/shared/scripts/Spec_OC_tot.nxs
<input type="checkbox"/>	18	corelli	vanadium_SA_file	/SNS/CORELLI/IPTS-21454/shared/scripts/SA_CC_tot.nxs
<input type="checkbox"/>	17	corelli	mask	[]
<input type="checkbox"/>	16	corelli	ub_matrix_file	/SNS/CORELLI/IPTS-21454/shared/250mK/optub_250mK.mat
<input type="checkbox"/>	15	arcs	e_max	0.95
<input type="checkbox"/>	14	arcs	e_step	0.01
<input type="checkbox"/>	13	arcs	e_min	-0.95
<input type="checkbox"/>	12	arcs	processed_vanadium	/SNS/ARCS/shared/autoreduce/vanadium_files/van155050.nxs
<input type="checkbox"/>	11	arcs	raw_vanadium	/SNS/ARCS/IPTS-24780/nexus/ARCS_155050.nxs.h5
<input type="checkbox"/>	10	arcs	mask	MaskBTPParameters.append({'Pixel': '1-7,122-128'}) MaskBTPParameters.append({'Pixel': '1-12,117-128', 'Bank': '70'}) MaskBTPParameters.
<input type="checkbox"/>	9	arcs	grouping	/SNS/ARCS/shared/autoreduce/ARCS_2X1_grouping.xml

CNCS Configuration

[home](#) > [cncs](#) > configuration

Configuring the automated reduction

Instrument team members can use this page to generate a new automated reduction script.

- Click the submit button to create a new automated reduction script.
- Click the reset to populate the form with default values.
- The `reduce_CNCS.py` will automatically be overwritten once you click the submit button.

List of parameters for CNCS reduction template:

Raw vanadium	<input type="text" value="/SNS/CNCS/IPTS-22728/nexus/CNCS_326713.nxs.h5"/>		
Processed vanadium	<input type="text" value="van_326713.nxs"/>		
Output directory	<input type="text"/>		
Vanadium integration	min <input type="text" value="49500.0"/>	max <input type="text" value="50500.0"/>	
Motor names	<input type="text" value="omega"/>		
Temperature names	<input type="text" value="SampleTemp,sampletemp,SensorB,SensorA,temp5,temp8,sensor0normal,Se"/>		
Grouping file	<input type="text" value="8 x 1"/>		
Create elastic nxspe	<input checked="" type="checkbox"/>		
Create MD nxsp	<input type="checkbox"/>		
Energy in meV	<input type="checkbox"/>		
Energy binning	E _{min} <input type="text" value="-0.1"/>	E _{step} <input type="text" value="0.005"/>	E _{max} <input type="text" value="0.95"/>
TOF offset	t ₀ <input type="text"/>	Auto-fit t ₀ to get E=0 at elastic peak <input type="checkbox"/>	
Time independent bck	min <input type="text"/>	max <input type="text"/>	
UB matrix	a <input type="text" value="1.0"/>	b <input type="text" value="1.0"/>	c <input type="text" value="1.0"/>
	alpha <input type="text" value="1.0"/>	beta <input type="text" value="1.0"/>	gamma <input type="text" value="1.0"/>
	u_vector <input type="text" value="1,0,0"/>	v_vector <input type="text" value="0,1,0"/>	
Masked Bank	Masked Tube	Masked Pixel	<input type="button" value="+"/>
<input type="text"/>	<input type="text"/>	<input type="text" value="121-128"/>	<input type="button" value="x"/>
<input type="text"/>	<input type="text"/>	<input type="text" value="1-8"/>	<input type="button" value="x"/>
<input type="text" value="36-50"/>	<input type="text"/>	<input type="text"/>	<input type="button" value="x"/>

Latest post-processing log entries for CNCS:

No recent changes

Example form

Uses standard Django fields →

```
class ReductionConfigurationCNCSForm(BaseReductionConfigurationForm):
    """
    Generic form for DGS reduction instruments
    """
    mask = forms.CharField(required=False, initial='')
    sub_directory = forms.CharField(required=False, initial='', widget=forms.TextInput(attrs={'class': 'font_resize'}))
    raw_vanadium = forms.CharField(required=False, initial='', widget=forms.TextInput(attrs={'class': 'font_resize'}))
    processed_vanadium = forms.CharField(required=False, initial='', widget=forms.TextInput(attrs={'class': 'font_resize'}))
    vanadium_integration_min = forms.FloatField(required=True, initial=84000)
    vanadium_integration_max = forms.FloatField(required=True, initial=94000)
    grouping = forms.ChoiceField(choices=[])
    e_pars_in_mev = forms.BooleanField(required=False)
    e_min = forms.FloatField(required=True, initial=-0.2)
    e_step = forms.FloatField(required=True, initial=0.015)
    e_max = forms.FloatField(required=True, initial=0.95)
    tib_min = forms.CharField(required=False, initial="", validators=[validate_float_list])
    tib_max = forms.CharField(required=False, initial="", validators=[validate_float_list])
    do_tib = forms.BooleanField(required=True)
    t0 = forms.CharField(required=False, initial="", validators=[validate_float_list])
    motor_names = forms.CharField(required=False, initial='huber,SERotator2,0xDilRot,CCR13VRot,SEOCRot,CCR10G2Rot,0x2WeldRot,ThreeSampleRot')
    temperature_names = forms.CharField(required=False, initial='SampleTemp,sampletemp,SensorC,SensorB,SensorA,temp5,temp8')
    create_elastic_nxspe = forms.BooleanField(required=False)
    create_md_nxs = forms.BooleanField(required=False)
    a = forms.FloatField(required=True, initial=7.76)
    b = forms.FloatField(required=True, initial=7.76)
    c = forms.FloatField(required=True, initial=7.02)
    alpha = forms.FloatField(required=True, initial=90)
    beta = forms.FloatField(required=True, initial=90)
    gamma = forms.FloatField(required=True, initial=90)
    u_vector = forms.CharField(required=False, initial="1,0,0", validators=[validate_float_list])
    v_vector = forms.CharField(required=False, initial="0,0,1", validators=[validate_float_list])
    auto_tzero_flag = forms.BooleanField(required=False)
```

List of fields used in the template →

```
# List of field that are used in the template
_template_list = ['mask', 'sub_directory', 'raw_vanadium', 'processed_vanadium', 'grouping',
                  'vanadium_integration_min', 'vanadium_integration_max',
                  'tib_min', 'tib_max', 'do_tib', 't0', 'motor_names', 'temperature_names',
                  'create_elastic_nxspe', 'create_md_nxs',
                  'alpha', 'beta', 'gamma', 'auto_tzero_flag',
                  'u_vector', 'v_vector', 'e_pars_in_mev',
                  'e_min', 'e_step', 'e_max', 'a', 'b', 'c']
```

If you need to populate drop-downs according to instrument, or initialize fields dynamically, do it here →

```
def __init__(self, *args, **kwargs):
    super(ReductionConfigurationCNCSForm, self).__init__(*args, **kwargs)

def set_instrument(self, instrument):
    """
    Populate instrument-specific options.
    @param instrument: instrument short name
    """
    self.fields['grouping'].choices = _get_choices(instrument)
```


The AMQ message

- The code executing the job is part of the post_processing_agent:

postprocessing/reduction_script_writer.py

- The web monitor doesn't have access to the file system, but the AR service does...

data_workflow/reporting/reduction/view_utils.py

```
def send_template_request(instrument_id, template_dict, user='unknown'):
    """
    Send an ActiveMQ message to request a new script

    @param instrument_id: Instrument object
    @param template_dict: dictionary of perproperties
    @param user: user that created the change
    """
    use_default = False
    if 'use_default' in template_dict:
        if type(template_dict['use_default']) == bool:
            use_default = template_dict['use_default']
        else:
            use_default = template_dict['use_default'].lower()=='true'

    encoded_dict = {}
    for key, value in template_dict.items():
        if isinstance(value, basestring):
            encoded_dict[key] = urllib.quote_plus(value)
        else:
            encoded_dict[key] = value

    # Send ActiveMQ request
    dasmon.view_util.add_status_entry(instrument_id,
                                     settings.SYSTEM_STATUS_PREFIX+'postprocessing',
                                     "Script requested by %s" % user)

    data_dict = {"instrument": str(instrument_id).upper(),
                 "use_default": use_default,
                 "template_data": encoded_dict,
                 "information": "Requested by %s" % user}
    data = json.dumps(data_dict)
    reporting_app.view_util.send_activemq_message(settings.REDUCTION_SCRIPT_CREATION_QUEUE, data)
    logging.info("Reduction script requested: %s", str(data))
```

The template

- Template resides in the instrument's shared area
- It can be modified by the instrument team

Replace keys
with values

```
#!/usr/bin/env python

#imports section
import sys, os, glob, filecmp, datetime, shutil
try:
    import ConfigParser
except:
    import configparser as ConfigParser
sys.path.append("/SNS/CNCS/shared/autoreduce")
from ARLibrary import * #note that ARLibrary would set mantidpath as well
sys.path.append("/opt/Mantid/bin")
from mantid.simpleapi import *
import numpy as np
import scipy.optimize as opt
import scipy.interpolate as interp

#parameters section
#this part changes with web input
MaskBTPParameters=[]
${mask}
#MaskBTPParameters.append({'Pixel': '1-43,95-128'})
#MaskBTPParameters.append({'Pixel': '1-7,122-128'})
#MaskBTPParameters.append({'Bank': '36-50'})#8T magnet
raw_vanadium="${raw_vanadium}"
processed_vanadium="${processed_vanadium}"
VanadiumIntegrationRange=[${vanadium_integration_min},${vanadium_integration_max}]#integration
grouping="${grouping}" #allowed values 1x1, 2x1, 4x1, 8x1, 8x2 powder
Emin="${e_min}"
Emax="${e_max}"
Estep="${e_step}"
E_pars_in_mev=${e_pars_in_mev}
TIB_min="${tib_min}"
TIB_max="${tib_max}"
doTIB=${do_tib}
T0="${t0}"
Motor_names="${motor_names}"
Temperature_names="${temperature_names}"
create_elastic_nxspe=${create_elastic_nxspe} #+-0.1Ei, 5 steps
create_MDnxs=${create_md_nxs}
a="${a}"
b="${b}"
c="${c}"
alpha="${alpha}"
beta="${beta}"
gamma="${gamma}"
uVector="${u_vector}"
vVector="${v_vector}"
sub_directory="${sub_directory}"
auto_tzero_flag = ${auto_tzero_flag}
```